

* In cattle, *Mycobacterium paratuberculosis* (Johne's disease) ELISA test results are interpreted as negative if less than or equal to 0.450, suspect if greater than 0.450 but less than 0.550, and positive if greater than or equal to 0.550. At the 0.550 cutoff, test sensitivity and specificity are estimated at 68% and 99%, respectively, in cattle.

In sheep and goats, *Mycobacterium paratuberculosis* (Johne's disease) ELISA test results are interpreted as negative if less than or equal to 0.100, suspect if greater than 0.100 but less than 0.300, and positive if greater than or equal to 0.300. At the 0.300 cutoff, test sensitivity and specificity are estimated at 51% and 95%, respectively, in small ruminants.

Brucella (B. melitensis) 3% Rose Bengal antibody test

Animal/Source	Specimen	Specimen Type	Results
1	(1) Melody	Serum	Negative
1	(2) Pepe	Serum	Negative
1	(3) Darlin	Serum	Negative
1	(4) Lady	Serum	Negative
1	(5) Rose	Serum	Negative
1	(6) Tessa	Serum	Negative
1	(7) Higgs	Serum	Negative
1	(8) Hazel	Serum	Negative
1	(9) Time	Serum	Negative
1	(10) Delilah	Serum	Negative
1	(11) Ginger	Serum	Negative
1	(12) Toby	Serum	Negative
1	(13) Badger	Serum	Negative
1	(14) Claire	Serum	Negative

Caprine arthritis encephalitis / ovine progressive pneumonia antibody cELISA

Animal/Source	Specimen	Specimen Type	Results
1	(1) Melody	Serum	Neg
1	(2) Pepe	Serum	Neg
1	(3) Darlin	Serum	Neg
1	(4) Lady	Serum	Neg
1	(5) Rose	Serum	Neg
1	(6) Tessa	Serum	Neg
1	(7) Higgs	Serum	Neg
1	(8) Hazel	Serum	Neg
1	(9) Time	Serum	Neg
1	(10) Delilah	Serum	Neg
1	(11) Ginger	Serum	Neg
1	(12) Toby	Serum	Neg
1	(13) Badger	Serum	Neg
1	(14) Claire	Serum	Neg

Corynebacterium pseudotuberculosis antibody - SHI

Animal/Source	Specimen	Specimen Type	Results
----------------------	-----------------	----------------------	----------------

1	(1) Melody	Serum	Neg @1:8
1	(2) Pepe	Serum	Neg @1:8
1	(3) Darlin	Serum	Neg @1:8
1	(4) Lady	Serum	Neg @1:8
1	(5) Rose	Serum	Neg @1:8
1	(6) Tessa	Serum	Neg @1:8
1	(7) Higgs	Serum	Neg @1:8
1	(8) Hazel	Serum	Neg @1:8
1	(9) Time	Serum	Neg @1:8
1	(10) Delilah	Serum	Neg @1:8
1	(11) Ginger	Serum	Neg @1:8
1	(12) Toby	Serum	Neg @1:8
1	(13) Badger	Serum	Neg @1:8
1	(14) Claire	Serum	Neg @1:8

M. paratuberculosis (Johnes Disease) Antibody ELISA

Animal/Source	Specimen	Specimen Type	Results
1	(1) Melody	Serum	0.00
1	(2) Pepe	Serum	0.027
1	(3) Darlin	Serum	0.001
1	(4) Lady	Serum	0.004
1	(5) Rose	Serum	0.016
1	(6) Tessa	Serum	0.00
1	(7) Higgs	Serum	0.001
1	(8) Hazel	Serum	0.00
1	(9) Time	Serum	0.001
1	(10) Delilah	Serum	0.00
1	(11) Ginger	Serum	0.005
1	(12) Toby	Serum	0.00
1	(13) Badger	Serum	0.00
1	(14) Claire	Serum	0.00